

UEV3 LDH domain

PDB:2I6T

Revision

Revision Type:created

Revised by:created

Revision Date:created

Entry Clone Accession:

Entry Clone Source:Open Biosystems

SGC Clone Accession: ubc81.171.471; plate SDC092:H7

Tag:N-terminal His-tag with integrated TEV-cleavage site: MHHHHHHSSGRENLYFQG.

Host:E.coli BL21 (DE3)

Construct

Prelude:

Sequence:

MHHHHHHSSGRENLYFQGSKSWANHENKTVNKITVVGGEGLGIACTLAISAKGIADRLVLLDLSEGTKGATMDLEIFNLPNVEISKD
LSASAHSKVVIPTVNSLGSSQSYLDVVQSNVDMFRALVPALGHYSQHSVLLVASQPVEIMTYVTWKLSTFPANRVIGIGCNLDSQRL
QYIITNVLKAQTSKGKEVWVIGEKGEDKVL TWSGQEEVVSHTSQVLSNRMELLRVKGQRSWSVGLSVADMVDSIVNNKKKVHVSVAL
AKGYYDINSEVFLSLPCILGTNGVSEVIKTTLKEDTVTEKLQSSASSIHS LQQQLKL

Vector:p28a-mhl

Growth

Medium:

Antibiotics:

Procedure:UEV3 LDH domain was expressed in E. coli BL21 (DE3) grown in Terrific Broth (Sigma) in the presence of 50 µg/ml of kanamycin at 37°C to an OD600 of 5.5-6. Protein expression was then induced by isopropyl-1-thio-D-galactopyranoside at a final concentration of 0.05 mM, and incubation continued overnight at 15°C. The culture was centrifuged; the cell pellets were collected, frozen with liquid nitrogen and stored at -80°C.

Purification

Procedure

Column 1: TALON metal-affinity resin column (BD Biosciences) .

Column 2: HighLoad 16/60 Superdex 200 column (GE Healthcare, Amersham)

The cleared lysate was loaded onto a TALON metal-affinity resin (BD Biosciences) column at 4°C. The column was washed with wash buffer A, wash buffer B and again wash buffer A, and the protein was eluted with elution buffer. The protein was further purified by gel filtration on a HighLoad 16/60 Superdex 200 column (GE Healthcare, Amersham) equilibrated with 20 mM Tris-HCl, pH 8.0, 0.5 M NaCl, 5% glycerol, 2 mM dithiothreitol and concentrated by ultrafiltration.

Extraction

Procedure

The cell pellet was resuspended in lysis buffer and lysed using Microfluidizer (peak pressure 18,000 psi). The lysate was cleared by centrifugation.

Concentration:

Ligand

MassSpec:

Crystallization: Purified UEV LDH domain was crystallized using the hanging-drop vapor-diffusion method. Crystals grew when the protein (6 mg/ml in the above gel-filtration buffer) was mixed with the reservoir solution in a 1:1 volume ratio, and the drop was equilibrated against a reservoir solution containing 26% PEG3350, 0.2 M ammonium sulfate, 0.1 M sodium cacodylate, pH 5.2 at 293K temperature.

NMR Spectroscopy:

Data Collection:

Data Processing: