

EIF5A

PDB:3CPF

Revision

Revision Type:created

Revised by:created

Revision Date:created

Entry Clone Accession:BC000751

Entry Clone Source:MGC AU5-C4

SGC Clone Accession:HPC076-B07

Tag:mhhhhhssgrenlyfq*g

Host:BL21-CodonPlus(DE3)-RIL

Construct

Prelude:

Sequence:

gSATFPMQCSALRKNGFVVLKGRPCKIVEMSTSKTGKHGHAKVHLVGIDIFTGKKYEDICPSTHNMDVPNIKRNDFQLIGIQDGYLS
LLQDSGEVREDLRLPEGDLGKEIEQKYDCGEEILITVLSAMTEEA AVA IKA

Vector:pET28-mhl (GI:134105571)

Growth

Medium:Terrific Broth

Antibiotics:

Procedure:LEX Bubbling. The target protein was expressed in E. coli by inoculating 100 mL of overnight culture grown in Luria-Bertani medium into a 1.8 L of Terrific Broth medium in the presence of 50 µg/mL kanamycin and 25 µg/mL chloramphenicol at 37 °C. When OD₆₀₀ reached ~3.0, the temperature of the medium was lowered to 15 °C and the culture was induced with 0.5 mM IPTG. The cells were allowed to grow overnight before harvested and flash frozen in liquid nitrogen and stored at -80 °C.

Purification

Procedure

The lysate was centrifuged at 15,000 rpm for 45 minutes and the supernatants were mixed with 5 mL 50% slurry of Ni-NTA beads and incubated at 4°C on rotary shaker for one hour. The mixture was then centrifuged at 2300 rpm for 5 min and the supernatant discarded. The beads were then washed with washing buffer containing 30 mM and 75 mM Imidazole, and finally the elution buffer. The flow-through was collected and further purified by a Superdex-75 gel filtration column pre-equilibrated with gel filtration buffer. Fractions were collected and digested with TEV protease. TEV protease was removed from the treated protein sample by adding 100 uL

50% slurry of Ni-NTA beads and the sample was purified with Superdex-75 gel filtration again. Fractions containing the protein were collected and concentrated with Amicon Ultra-15 centrifugal filter. The purity of the preparation is tested by SDS-PAGE to be around 99%.

Extraction

Procedure

Frozen cells from 1.8L TB culture were thawed and resuspended in 150 mL extraction buffer with freshly added 0.5% CHAPS and 2mM BME, and supplemented with protease inhibitor cocktail (SIGMA Catalog # P8849), and 3 μ L benzonase (Sigma Catalog # E1014, 250U/ μ L), and lysed using microfluidizer at 15,000 PSI.

Concentration: 29.5 mg/mL

Ligand

N/A **MassSpec:** Native: 15149.67, expected 15149.41

Crystallization: Crystallization was setup using sitting drops with Red Wings and SGC-I screens initially. Small crystals were seen at condition SC08.

Optimization was done using hanging drop vaporization.

Crystal used for structure determination were grown in: 22.0% PEG3350, 0.2M $(\text{NH}_4)_2\text{SO}_4$, 0.1 M NaCac pH 5.5. Cryo used 20% PEG 3350 + 20% EG.

Last updated by ytong 20080410

NMR Spectroscopy:

Data Collection:

Data Processing: