

RAB3B

PDB:3DZ8

Revision

Revision Type:created

Revised by:created

Revision Date:created

Entry Clone Accession:BC005035

Entry Clone Source:MGC AU35-H11

SGC Clone Accession:HPC085-A01

Tag:mhhhhhssgrenlyfq*g

Host:BL21-CodonPlus(DE3)-V2R

Construct

Prelude:

Sequence:

mhhhhhssgrenlyfqgNFDYMFKLLIIGNSSVGKTSFLFRYADDTFTPAFVSTVGIDFKVKTVYRHEKRVKLQIWDTAGQERYRT
ITTAYYRGAMGFILMYDITNEESFNAVQDWATQIKTYSWDNAQVILVGNKCDMEEERVVPTEKGQLLAEQLGFDFFEASAKENISVR
QAFERLVDAICDKMSDS

Vector:pET28-mhl (GI:134105571)

Growth

Medium:Terrific Broth

Antibiotics:

Procedure:LEX Bubbling. The target protein was expressed in E. coli by inoculating 100 mL of overnight culture grown in Luria-Bertani medium into a 1.8 L of Terrific Broth medium in the presence of 50 µg/mL kanamycin and 25 µg/mL chloramphenicol at 37 °C. When OD₆₀₀ reached ~3.0, the temperature of the medium was lowered to 15 °C and the culture was induced with 0.5 mM IPTG. The cells were allowed to grow overnight before they were harvested and flash frozen in liquid nitrogen and stored at -80 °C.

Purification

Procedure

The lysate was centrifuged at 16,000 rpm for 60 minutes and the supernatant was mixed with 3 mL 50% Ni-NTA beads, and incubated at 4 °C for 2 hours on a roller drum. The supernatant was then passed through a gravity column (Poly-Prep, Bio-Rad, Catalog #731-1550) and the beads were washed using 50 mL washing buffer twice. The protein bound to beads were eluted using 15 mL elution buffer. The flow-through was collected and loaded onto Superdex-75 26/60 gel filtration column. Eluted fractions were pooled and concentrated using amicon centrifugal

filter (m.w. cut-off 10,000). The purity of the proteins was higher than 95% judged by SDS-PAGE.

Extraction

Procedure

Frozen cells from 3.6L TB culture were thawed and resuspended in 400 mL extraction buffer with freshly added 0.5% CHAPS and 10U/mL Benzonase (Sigma Catalog # E1014, 250U/ μ L), and supplemented with 1mL protease inhibitor cocktail (SIGMA Catalog # P8849), and lysed using sonication at 10 seconds 50% duty cycle for 8 minutes at 120W.

Concentration: 18.0 mg/mL

Ligand

GDPMassSpec: Native: 22155.29, expected 22154.94

Crystallization: Crystallization was setup using sitting drops with Red Wings and SGC-I screens initially. 2D-plate xtals were seen at condition SG04, SD11 several days after setup.

Crystal used for data collection was grown at 2.0M (NH₄)₂SO₄, 0.2M NaAc, 0.1M HEPES pH 7.5, 5% MPD.

Last updated by y tong 20080730

NMR Spectroscopy:

Data Collection:

Data Processing: