

# CYP51A1

PDB:3I3K

## Revision

**Revision Type:**created

**Revised by:**created

**Revision Date:**created

**Entry Clone Accession:**AT52-A8 (BC032322)

**Entry Clone Source:**MGC

**SGC Clone Accession:**

**Tag:**N-terminal: MAKKT;

C-terminal: 6His-tag

**Host:***E. coli* JM109 (Stratagene).

## Construct

**Prelude:**

**Sequence:**

LPAGVKSPPIYIFSPIPLGHAIAFGKSPIEFLENAYEKYGPVFSFTMVGKTFTYLLGSDAAALLFNSKNEDLNAEDVYSRLTTPVFG  
KGVAYDVPNPVFLQKKMLKSGLNIAHFKQHVSIIKETKEYFESWGESGEKNVFEALSELIILTASHCLHGKEIRSQLNEKVAQLY  
ADLDGGFSHAAWLLPGWLPLPSFRRRDRAHREIKDIFYKAIQKRRQSQEKIDDILQTL LDATYK DGRPLTDDEVAGMLIGLLLAGQH  
TSSTTSAWMGFFLARDKTLQKKCYLEQKTVCGENLPPLTYDQLKDLNLLDRCIKETLRLRPPVMIMMRMARTPQTVAGYTIPP GHQV  
CVSPTVNQRLKDSWVERLDFNPDRYLQDNPASGEKFAYVPFGAGRHR CIGENFAYVQIKTIWSTMLRLYEFDLIDGYFPTVNYTTMI  
HTPENPVIRYKRR

**Vector:**pCW-LIC-29

## Growth

**Medium:**

**Antibiotics:**

**Procedure:**CYP51A1 was co-expressed with GroEL/ES in *E. coli* JM109 in TB medium. Cells were grown at 37 degC to an OD600 of 1.0 and induced by 0.5mM IPTG and 4mg/ml of arabinose and in the presence of 0.5mM  $\delta$ -aminolevulinic acid and incubated 48 hours at 26°C.

## Purification

### Procedure

The lysate was centrifuged at 60.000g for 60min. The supernatant was loaded onto 5ml NiHiTrap column (Amersham Biosciences) equilibrated with buffer A. The column was washed with buffer A and protein was eluted using a linear gradient of 5-100% Buffer B. The protein was further purified by ion-exchange chromatography on SourceS column (Amersham Biosciences), equilibrated with buffer 5mM KPi, pH 7.4, 20% glycerol and eluted with linear gradient of Buffer C.

## Extraction

### Procedure

Collected/resuspended cells (50mM potassium phosphate, 300mM NaCl, 20% glycerol, pH 7.4, 0.4mM PMSF) were disrupted in a high-pressure Microfluidizer (Microfluidics Corp.) at 18.000 psi.

**Concentration:** 20 mg/ml

### Ligand

**MassSpec:** Expected MW is 52550, measured mass is 52551.

**Crystallization:** Purified CYP51A1 was crystallized in presence of ketoconazole using hanging drop vapor diffusion method drop at 18 °C by mixing 1µl of the protein solution with 1µl of the reservoir solution containing 0.1 M Hepes pH 7.5, 2.5 M Ammonium Sulfate.

### NMR Spectroscopy:

#### Data Collection:

#### Data Processing: