

PKM2

PDB:3U2Z

Revision

Revision Type:created

Revised by:created

Revision Date:created

Entry Clone Accession:BC007952

Entry Clone Source:MGC: AU36-H1

SGC Clone Accession:HPC002-A01

Tag:N-terminal His6-tag, not removed

Host:BL21-V2R

Construct

Prelude:PANK3: M1-P531

Sequence:

mgsshhhhhssglvprgsMSKPHSEAGTAFIQTQQLHAAMADTFLEHMCRLDIDSPITARNTGIICTIGPASRSVETLKEMIKSG
MNVARLNFSGHGTHEYHAETIKNVRTATESFASDPILYRPVAVALDTKGPEIRTGLIKSGTAEVELKKGATLKITLDNAYMEKCDEN
ILWLDYKNICKVVEVGSKIYVDDGLISLQVKQKGADFLVTEVENGGSLGSKKGVNLPGAAVDLPVSEKDIQDLKFGVEQDVMVFA
SFIRKASDVHEVRKVLGEKGKNIKIISKIENHEGVRRFDEILEASDGMVARGDLGIEIPA EKVF LAQKMMIGRCNRAGKPVICATQ
MLESMIKKPRPTRAEGSDVANAVLDGADCIMLSGETAKGDYPLEAVRMQH LIAREAEAAIYHLQLFEELRR LAPITSDPTEATAVGA
VEASFKCCSGAIIVLTKSGRSAHQVARYRPRAPIIAVTRNPQTARQAHL YRGIFPVLCKDPVQEAWAEDVDLRVNFAMNVGKARGFF
KKGDVVIVLTGWRPGSGFTNTMRVVPVP

Vector:pET28a-LIC

Growth

Medium:

Antibiotics:

Procedure:The seeds were grown in 80 mL Luria-Bertani broth media supplemented with 50 µg/mL kanamycin at 37 °C overnight. The following morning, all of the seeds were inoculated 1800 mL of Terrific Broth media supplemented with 50 µg/mL kanamycin, 8 g/l glycerol and approximately 500 µl antifoam in glass flasks in the Large Scale Expression System (LEX). Cells were grown at 37 °C until OD_{600nm} of 4.0 and were then induced by addition of IPTG to a concentration of 0.5 mM. Protein expression was allowed to continue over night at 18 °C.

Purification

Procedure

The supernatant was passed through DE52 (Whatman) column equilibrated with the binding buffer and then loaded onto 3 mL Ni-NTA column (Qiagen) equilibrated with the same binding buffer at 4 °C. The Ni-NTA column was washed with 150 mL of the wash buffer and the protein

was eluted with 15 mL of the elution buffer. The eluate was further purified by size-exclusion chromatography (Superdex 200) equilibrated with 10 mM HEPES pH 7.5, 150mM KCl, 2mM TCEP, 5% glycerol, and 5mM MgCl₂. The protein was concentrated using an Amicon Ultra centrifugal filter to the concentration of 50 mg/mL.

Extraction

Procedure

Cultures were centrifuged and the cell pellets were suspended in 100 mL of the binding buffer with a protease inhibitor cocktail (0.1 mM M benzamidine-HCl and 0.1 mM phenylmethyl sulfonyl fluoride) and flash frozen. The thawed cell pellet was lysed by a combination of 0.5% CHAPS (Sigma) and sonication. The lysate was centrifuged at 15000 rpm for 30 min and the supernatant was used for subsequent steps of purification.

Concentration: 50 mg/mL

Ligand

6-(3-aminobenzyl)-4-methyl-2-methylsulfinyl- 4,6-dihydro-5H-thieno[2',3':4,5]pyrrolo[2,3-d]pyridazin-5-one

MassSpec:

Crystallization: Crystallization trials were set up using the vapor diffusion method and the protein drop was equilibrated against a reservoir solution with 1:1 volume ratio. Prior to crystallization, the purified protein was incubated overnight at 4 °C in the presence of the activator (to 5 ~10 mM final concentration). Crystals of activator bound PKM2 were grown at 25% PEG-3350, 0.1M ammonium sulfate, 0.1M Bis-Tris, pH6.5.

NMR Spectroscopy:

Data Collection:

Data Processing: