

PRMT6

PDB:4HC4

Revision

Revision Type:created

Revised by:created

Revision Date:created

Entry Clone Accession:GI:227908866

Entry Clone Source:MGC

SGC Clone Accession:HRMT1L6:PBC001-A12:C213642

Tag:N-terminal: His-tag with integrated TEV protease site:MHHHHHHSSGRENLVFQ*G

Host:Sf9

Construct

Prelude:

Sequence:

gMSQPKRRKLESGGGGEGGEGTEEDGAEREAALEPRRTKRERDQLYYECYSDVSVHEEMIADRVRTDAYRLGILRNWAALRGKTV
LDVGAGTGILSIFCAQAGARRVYAVEASAIWQQAREVVRFNGLEDRVHVLPGPVETVELPEQVDAIVSEWMGYGLLHESMLSSVLHA
RTKWLKEGGLLLPASAEFIAPISDQMLEWRLGFWSQVKQHYGVDMSCLEGFATRCLMGHSEIVVQGLSGEDVLARPQRFAQLELSR
AGLEQELEAGVGGRFRCSCYGSAPMHGFAIWFQVTFPGGESEKPLVLSTSPFHPATHWKQALLYLNPEVQVEQDQTDVSGEITLLPSR
DNPRLRLVLLRYKVGDQEEKTKDFAMED

Vector:pFBOH-MHL

Growth

Medium:HRMT1L6 was expressed in Sf9 cells

Antibiotics:

Procedure:

Purification

Procedure

For purification the cell paste was thawed and resuspended in lysis buffer containing 20 mM Tris-HCl, pH 8.0, 500 mM NaCl, 5 mM imidazole, 2 mM β -mercaptoethanol, 5% glycerol, 0.6% NP-40, protease inhibitor cocktail (Roche), 3000U of benzonase (Novagen). Cells were lysed by brief sonication. The clarified lysate was loaded onto a 2 mL TALON column (Clontech). The column was washed with 50 column volumes of 20 mM Tris-HCl buffer, pH 8.0, containing 500 mM NaCl, 5% glycerol and 5 mM imidazole, and the protein was eluted with elution buffer (20 mM Tris-HCl, pH 8.0, 250 mM NaCl, 5% glycerol, 250 mM imidazole). The eluted protein was loaded on a Superdex200 column (GE Healthcare), equilibrated with 20 mM Tris-HCl buffer, pH 8.0, and 150 mM NaCl. Pooled fractions containing HRMT1L6 were subjected to TEV treatment

to remove His-tag. The protein was further purified to homogeneity by ion-exchange chromatography.

Enzymatic treatment: TEV

Extraction

Procedure

Cells were harvested by centrifugation at 5,000rpm. The cell pellets were frozen in liquid nitrogen and stored at -80°C.

Concentration: 16 mg/ml

Ligand

MassSpec: expected mass is 42022.6 Da, measured mass is 42132.48Da.

Crystallization: Purified HRMT1L6 (4.9 mg/ml) was complexed with (SAH, Sigma) at 1:5 molar ratio of protein:SAH and crystallized using the hanging drop vapor diffusion method at 20 °C by mixing 1 µl of the protein solution with 1 µl of the reservoir solution containing 15% PEG3350, 0.1 M succinate acid, pH 7.0.

NMR Spectroscopy:

Data Collection:

Data Processing: