

**Target ID:**MYO5BA

**Entry Clone ID:** MYO5BA-s001

**Allele ID:** MYO5BA-a002

**Construct ID:** MYO5BA-c002

**Clone ID:** MYO5BA-k002

**Expression ID:** MYO5BA-e007

**Purification ID:** MYO5BA-p003

**Entry clone source:** MGCA

**Vector:** Brazil-1-pNIC28-Bsa4

**E.coli strain:** BL21(DE3)-R3-pRARE2

**Tags and additions:** N-terminal, TEV protease cleavable hexahistidine tag, tag not been cleaved

### Coding DNA sequence

```
ATGAAGAAAGCCCAGGACCTAGAAGC
TGCCCAGGCATTGCCCCAGAGTGAGA
GGAAGCGCCATGAGCTCAACAGGCAG
GTCACGGTCCAGCGGAAAGAGAAGGA
TTTCCAGGGCATGCTGGAGTACCACA
AAGAGGACGAGGCCCTCCTCATCCGG
AACCTGGTGACAGACTTGAAGCCCCA
GATGCTGTCGGGCACAGTGCCCTGTC
TCCCCGCCTACATCCTCTACATGTGC
ATCCGGCACGCGGACTACACCAACGA
CGATCTCAAGGTGCACTCCCTGCTGA
CCTCCACCATCAACGGCATTAAGAAA
GTCCTGAAAAAGCACAATGATGACTT
TGAGATGACGTCATTCTGGTTATCCA
ACACCTGCCGCCTTCTTCACTGTCTG
AAGCAGTACAGCGGGGATGAGGGCTT
CATGACTCAGAACTGCAAAGCAGA
ATGAACACTGTCTTAAGAATTTTGAC
CTCACCGAATACCGTCAGGTGCTGAG
TGACCTTTCCATTCAGATCTACCAGC
AGCTCATTAATAATTGCCGAGGGCGTG
TTACAGCCGATGATAGTTTCTGCCAT
GTTGGAAAATGAGAGCATTCAAGGTC
TATCTGGTGTGAAGCCCACCGGCTAC
CGGAAGCGCTCCTCCAGCATGGCAGA
TGGGGATAACTCATACTGCCTGGAAG
CTATCATCCGCCAGATGAATGCCTTT
CATAAGTCATGTGTGACCAGGGCTT
GGACCCTGAGATCATCCTGCAGGTAT
TCAAACAGCTCTTCTACATGATCAAC
GCAGTGAATCTTAACAACCTGCTCTT
GCGGAAGGACGTCTGCTCTTGAGCA
CAGGCATGCAACTCAGGTACAATATA
AGTCAGCTTGAGGAGTGGCTTCGGGG
AAGAAACCTTCACCAGAGTGGAGCAG
TTCAGACCATGGAACCTCTGATCCAA
GCAGCCCAGCTCCTGCAATTAAAGAA
```

GAAAACCCAGGAGGACGCAGAGGCTA  
TCTGCTCCCTGTGTACCTCCCTCAGC  
ACCCAGCAGATTGTCAAAATTTTAAA  
CCTTTATACTCCCCTGAATGAATTTG  
AAGAACGGGTAACAGTGGCCTTTATA  
CGAACAATCCAGGCACAACACTACAAGA  
GCGGAATGACCCTCAGCAACTGCTAT  
TAGATGCCAAGCACATGTTTCCTGTT  
TTGTTTCCATTTAATCCATCTTCTCT  
AACCATGGACTCAATCCACATCCCAG  
CGTGTCTCAATCTGGAATTCCTCAAT  
GAAGTCTGA

### **Final protein sequence**

mhhhhhssgvdlgtenlyfq\*smLN  
RQVTVQRKEKDFQGMLEYHKEDEALL  
IRNLVTDLKPQMLSGTVPCLPAYILY  
MCIRHADYTNDDLKVHSLTSTINGI  
KKVLKKHNDDFEMTSFWLSNTCRLH  
CLKQYSGDEGFMTQNTAKQNEHCLKN  
FDLTEYRQVLSDSLQIYQQLIKIAE  
GVLQPMIVSAMLENESIQLSGVKPT  
GYRKRSSSMADGDNSYCLEAIIRQMN  
AFHTVMCDQGLDPEIILQVFKQLFYM  
INAVTLNNLLLRKDVCSWSTGMQLRY  
NISQLEEWLRGRNLHQSGAVQTMEPL  
IQAAQLLQLKKKTQEDAEAICSLCTS  
LSTQQIVKILNLYTPLNEFEERVTV  
FIRTIQAQLQERNDPQQLLLDAKHM  
FVLFPPNPSSLTMDSIHIPACLNLEF  
LNEV

Sequence MHHHHHHSSGVDLG TENLYFQ\*SM is the purification tag (lower case) plus TEV protease recognition site (\*) .

### **Expression**

#### **Expression strain**

BL21(DE3)-R3-pRARE2

#### **Transformation**

The construct DNA was transformed into competent cells of the expression strain by a standard heat shock procedure.

#### **Glycerol stock preparation**

One colony from the transformation was used to inoculate 1 ml of LB media containing 50  $\mu$ g/ml kanamycin and 34  $\mu$ g/ml chloramphenicol, which was placed in a 37 $^{\circ}$ C shaker overnight. The next day glycerol stocks were prepared from this overnight culture.

### **Expression**

5 ml glycerol stock were used to inoculate 50 ml of LB media containing 50  $\mu\text{g/ml}$  kanamycin and 34  $\mu\text{g/ml}$  chloramphenicol, which was placed in a 37°C shaker overnight. The next day this starter culture was used to inoculate 6L of LB media (1 ml starter culture used per 1L) containing 50  $\mu\text{g/ml}$  kanamycin. When the OD<sub>600</sub> reached approximately 1.0 the temperature was reduced to 18°C and the cells were induced by the addition of 0.2 mM IPTG. The expression was continued overnight at 18°C.

### **Cell harvest**

Cells were harvested by centrifugation at 6500 rpm for 11 min at 4°C after which the supernatant was poured out and the cell pellet placed in a -20°C freezer.

### **Purification**

#### **Cell Lysis**

Cell pellets from 2 liter expression were slowly thawed on ice. Afterwards the cell pellets were dissolved in approximately 80 ml binding buffer and broken by sonication for 15 min with on/off settings of 10 sec and 35% amplitude. After lyses the solids were separated from the supernatant by centrifugation at 4°C for 60 min at 50,000 xg. The clear supernatant was transferred to a fresh 50 ml Falcon tube for further purification.

binding buffer: 25 mM Hepes pH 7.5, 500 mM NaCl, 5% Glycerol, 10 mM Imidazole, 0.5 mM TCEP, 1:500 Protease inhibitor cocktail Set III

#### **Column 1**

Ni-NTA (4 ml volume in a gravity-flow column).

The clarified cell extract was further purified on a 4 ml of Ni-NTA column. The supernatant, already containing binding buffer, was applied on the column twice before washing and eluting. During the washing step small 5 ml portions of washing buffer were added ten times consecutively. The protein was eluted with four times 5 ml of Elution Buffer.

Binding Buffer: see above

Wash Buffer: 20 mM Hepes pH 7.5, 500 mM NaCl, 5% Glycerol, 30 mM Imidazole

Elution Buffer: 20 mM Hepes pH 7.5, 500 mM NaCl, 5% Glycerol, 500 mM Imidazole

#### **Column2**

Superdex 200 10/300 column

The wash buffer fractions and elution buffer fractions from column 1 were pooled separately and concentrated to 5 ml with a 30 kDa mwco spin concentrator. The column had been pre-equilibrated with gel filtration buffer at 1.0 ml/min. The 5 ml protein sample was injected onto the column and 1.0 ml fractions were collected. The protein eluted at between 85 ml and 100 ml volume.

Gel Filtration Buffer: 10 mM HEPES, pH 7.5; 200 mM NaCl; 5% Glycerol; 0.5 mM TCEP

### **Concentration**

The eluted protein was concentrated to 10.7 mg/ml and stored at -80°C.

## Crystallization

Before crystallization protein was treated with Trypsin to allow mild tryptic digest (10mg/mL of Trypsin per 1mg/mL of protein). Crystals were grown by vapour diffusion in sitting drop at 4°C. by setting up 12.6 mg/ml of protein in the presence of 5 mM NADP<sup>+</sup>. Crystals appeared in a sitting drop consisting of 100 nl protein and 50 nl well solution which had been equilibrated against 20 ml well solution containing 0.2 M sodium nitrate and 20% (w/v) PEG 3350. Crystals were mounted in the presence of 25% ethylene glycol and flash cooled in liquid nitrogen.

## Data collection

Resolution: 2.25 Å

X-ray source: Diamond Light Source beamline IO2