

TBL1XR1

PDB:4LG9

Revision

Revision Type:created

Revised by:created

Revision Date:created

Entry Clone Accession:NP_078941.2.

Entry Clone Source:MGC 11-I9 (NP_078941.2.)

SGC Clone Accession:TBL1XR1_02b (APC048:A06): E134-K514

Tag:N-terminal tag: MGSSHHHHHHSSGLVPLGS

Host:baculovirus

Construct

Prelude:

Sequence:

MGSSHHHHHHSSGLVPRGSEENGAHTIANNHTDMMEVDGDVEIPPNAKAVVLRGHESEVFICAWNPNVSDLLASGSGDSTARIWNLSEN
STSGSTQLVLRHCIREGGQDVPSNKDVTSLDWNSEGTLLATGSYDGFARIWTKDGNLASTLGQHKGPWFALKWNKKGNFILSAGVDK
TTIIWDAHTGEAKQQFPFHSAPALVDWQSNNTFASCSTDMCIHVCKLGQDRPIKTFQGHTNEVNAIKWDP TGNLLASCSDDMTLKI
WSMKQDNCVHDLQAHNKEIYTIKWSPTGPGTNNPNANLMLASASFDSTVRLWDVDRGICHTMTKHQEPVYSVAFSPDGRYLASGSF
DKCVHIWNTQTGALVHSYRGTGGIFEVCWNAAGDKVGASASDGSVCVLDLRK

Vector:pFBOH-Lic

Growth

Medium:

Antibiotics:

Procedure:Baculovirus P1, P2, P3

Purification

Procedure

IMAC: Unclarified lysate was mixed with 1-2 mL of Ni-NTA superflow Resin (Qiagen) per 40 mL lysate. The mixture was incubated with mixing for at least 45 minutes at 4°C. The mixture was then loaded onto an empty comLum (BioRad) and washed with 100 mL wash buffer.

Samples were eluted from the resin by exposure to 2-3 column volumes (approx. 10-15 mL) of elution buffer. Concentration of eluted protein was estimated by OD280

Gel filtration chromatography: An XK 26x65 column (GE Healthcare) packed with HighLoad Superdex 75 resin (GE Healthcare) was pre-equilibrated with gel filtration buffer for 1.5 column volumes using an AKTA explorer (GE Healthcare) at a flow rate of 1.0 mL/min. The dialyzed sample from the IMAC step (approx. 15 mL) was loaded onto the column at 1.5 mL/min, and

2mL fractions were collected into 96-well plates (VWR 40002-012) using peak fractionation protocols). Fractions observed by a UV absorption chromatogram to contain the protein were pooled.

Extraction

Procedure

Frozen cell pellet contained in bags (Beckman 369256) obtained from 6L of culture were thawed by soaking in warm water. Each cell pellet was resuspended in 25-40 mL lysis buffer and homogenized using an Ultra-Turrax T8 homogenizer (IKA Works) at maximal setting for 30-60 seconds per pellet. Cell lysis was accomplished by sonication (Virtis408912, Virsonic) on ice: the sonication protocol was 10 sec pulse at half-maximal frequency (5.0), 10 second rest, for 10 minutes total sonication time per pellet.

Concentration: Purified proteins were concentrated using 15 mL concentrators with a 5,000 molecular weight cut-off (Amicon Ultra-15, UFC900524, Millipore) at 3750 rpm, typically resulting in a final concentration around 20 mg/mL.

Ligand

MassSpec:

Crystallization: Recombinant human TBL1XR1 WD40 repeats was crystallized using the hanging drop vapour diffusion method at 18 °C. The crystals were obtained in a buffer containing 0.1M Sodium Hepes, pH 7.5, 2.0 M Ammonium Sulphate, 2% PEG 400. Crystals were soaked in a cryoprotectant consisting of 100% reservoir solution and 12% glycerol.

NMR Spectroscopy:

Data Collection:

Data Processing: