

# PRMT6

**PDB:**5E8R

## Revision

**Revision Type:**created

**Revised by:**created

**Revision Date:**created

**Entry Clone Accession:**GI:227908866

**Entry Clone Source:**MGC

**SGC Clone Accession:**

**Tag:**N-terminal: His-tag with integrated TEV protease site: MHHHHHHSSGRENL<sup>Y</sup>FQ\*G

**Host:**Sf9

## Construct

**Prelude:**

**Sequence:**

gMSQPKRRKLESGGGGEGEGTEEEEDGAEREAAALERPRTKRERDQLYYECYSDVSVHEEMIADRVRTDAYRLGILRNWAALRGKTV  
LDVGAGTGILSIFCAQAGARRVYAVEASAIWQQAREVVRFNGLEDRVHVLPGPVETVELPEQVDAIVSEWMGYGLLHESMLSSVLHA  
RTKWLKEGGLLLPASAELFIAPISDQMLEWRLGFWSQVKQHYGVDMSCLEGFATRCLMGHSEIVVQGLSGEDVLARPQRFAQLELSR  
AGLEQELEAGVGGRFRCSYGSAPMHGFAIWFQVTFPGGESEKPLVLSTSPFHPATHWKQALLYLNPEVQVEQD<sup>T</sup>DVSGEITLLPSR  
DNPRLRLVLLRYKVGDQEEKTKDFAMED

**Vector:**pFBOH-MHL

## Growth

**Medium:**

**Antibiotics:**

**Procedure:**PRMT6 was expressed in Sf9 cells

## Purification

### Procedure

For purification the cell paste was thawed and resuspended in lysis buffer containing 20 mM Tris-HCl, pH 8.0, 500 mM NaCl, 5 mM imidazol, 2 mM  $\beta$ -mercaptoethanol, 5% glycerol, 0.6% NP-40, protease inhibitor cocktail (Roche), 3000 U of benzonase (Novagen). Cells were lyzed by brief sonication. The clarified lysate was loaded onto a 2 mL TALON column (Clontech). The column was washed with 50 column volumes of 20 mM Tris-HCl buffer, pH 8.0, containing 500 mM NaCl, 5% cglycerol and 5 mM imidazole, and the protein was eluted with elution buffer (20 mM Tris-HCl, pH 8.0, 250 mM NaCl, 5% glycerol, 250 mM imidazole). The eluted protein was loaded on a Superdex200 column (GE Healthcare), equilibrated with 20 mM Tris-HCl buffer, pH 8.0, and 150 mM NaCl. Pooled fractions containing HRMT1L6 were subjected to TEV treatment

to remove His-tag. The protein was further purified to homogeneity by ion-exchange chromatography.

## **Extraction**

### **Procedure**

Cells were harvested by centrifugation at 5,000 rpm. The cell pellets were frozen in liquid nitrogen and stored at -80°C.

**Concentration:** 19.6 mg/ml - Enzymatic treatment: TEV

### **Ligand**

**MassSpec:** expected mass is 42022.6 Da, measured mass is 42132.48 Da.

**Crystallization:** Purified HRMT1L6 (5.6 mg/ml) was mixed with S-adenosyl-L-homocysteine (SAH, Sigma) and MS023 at 1:5 molar ratio of protein:compound and crystallized using the sitting drop vapor diffusion method at 20 °C by mixing 1 µl of the protein solution with 1 µl of the reservoir solution containing 10% PEG 3350, 0.2 M MgCl<sub>2</sub>, 0.1 M sodium cacodylate, pH 5.6. The reservoir solution containing 10% ethylene glycol was used as cryo protectant.

### **NMR Spectroscopy:**

#### **Data Collection:**

#### **Data Processing:**