

# HMCES

**PDB:**5KO9

## Revision

**Revision Type:**created

**Revised by:**created

**Revision Date:**created

**Entry Clone Accession:**BC010125

**Entry Clone Source:**MGC AU62-F2

**SGC Clone Accession:**PBC006-B11

**Tag:**C-terminal His6-tag, removed

**Host:**BL21(DE3)V2R-pRARE2

## Construct

**Prelude:**

**Sequence:**

CGRTSCHLPRDVLTRACAYQDRRGQRLPEWRDPDKYCPSYNKSPQSNPVLRLHFEKDADSSERIAPMRWGLVPSWFKESDPS  
KLQFNTTNCRSDTVMEKRSEFKVPLGKGRRCVVLADGFYEWQRCQGTNQRQPYFIYFPQIKTEKSGSIGAADSPENWEKVWDNRLLT  
MAGIFDCWEPPEGDDVLYSYTITVDSCGLSDIHRMPAILDGEEAVSKWLDGFEVSTQEALKLIHPTENITFHAVSSVVNNSRNN  
TPECLAPVaenlyfq DNA sequence has been verified by sequencing

**Vector:**pNIC-CH

## Growth

**Medium:**

**Antibiotics:**

**Procedure:**LEX Bubbling - For native protein: The target protein was over-expressed in E. coli at 37degree C by inoculating 20 mL of overnight culture grown in Luria-Bertani medium into 2 L TerrificBroth medium in the presence of 40 ug/mL kanamycin and 25 ug/mL chloramphenicol. When the OD600 of the culture reached ~1.5, the temperature was lowered to 18 degree and the culture was induced with 0.25 mM final IPTG concentration. The cells were allowed to grow overnight before harvested by centrifugation (5,000 rpm Beckman JLA-8.1000rotor 15 min) and flash frozen in liquid nitrogen and stored at -80 degree.

## Purification

**Buffers**

Washing Buffer: 50 mM Tris-HCl pH8.0, 500 mM NaCl, 50 mM imidazole - Elution Buffer: 50 mM Tris-HCl pH8.0, 500 mM NaCl, 5% Glycerol, 250 mM imidazole - Dialysis Buffer: 50 mM Tris-HCl pH8.0, 300 mM NaCl, 2 mM TCEP - Ion Exchange Buffer A(B): 20 mM Tris-HCl pH 8.0, 0(1) M NaCl, 1 mM TCEP - Gel Filtration Buffer: 20 mM Tris-HCl pH 8.0, 150 mM NaCl, 1

mM TCEP

#### **Procedure**

The lysate was centrifuged at 15,000 rpm for 30 minutes. The supernatant was supplemented with 4 mL Ni-NTA resin (50% slurry) and incubated on a rotary drum for 1 hour at 4 degree, then loaded onto Bio-rad gravity column. The beads deposited in the open column was then washed with 20CV lysis buffer followed by 10CV washing buffer. Bound proteins were eluted using 5CV elution buffer. The N-terminal His-tag was removed by overnight incubation with TEV protease (1:30 w/w) at 4 degree during dialysis against the dialysis buffer. Uncut proteins and TEV protease were removed by passing the solution through 1mL Ni-NTA beads, and the target protein was further purified by anion-exchange chromatography on a 5mL HiTrap Q column (GE Healthcare). The protein was further purified using gel filtration on a HighLoad 16/60 Superdex 200 column (GE Healthcare) pre-equilibrated with gel filtration buffer. Fractions containing the target protein were pooled and concentrated by centrifugal filters (Amicon mwco 10kDa). The yields of the protein was about 4.2 mg per litre bacterial culture with >95% purity.

#### **Extraction**

##### **Buffers**

50 mM Tris-HCl pH 8.0, 500 mM NaCl, 5% Glycerol

##### **Procedure**

2L native cell pellet was resuspended in a total volume of 200 ml extraction buffer with 1mM PMSF/Benzamidine freshly added and the cells disrupted by sonication for 10 mins at 5" on 10" off duty cycle at 108W output power.

**Concentration:** Concentration used for crystallization : 18.0 mg/mL

##### **Ligand**

**MassSpec:** measured (after His-tag cleavage): 31.632.9 Da

**Crystallization:** Crystals used for structure determination were grown at 293K in sitting drop Vapor Diffusion by mixing 1 uL protein solution with 1 uL well solution containing 0.1M BTP, 2% Tacsimate, 20% PEG3350, pH 6.5. The crystals were cryoprotected by 30% Glycerol before flash frozen in liquid nitrogen.

##### **NMR Spectroscopy:**

##### **Data Collection:**

##### **Data Processing:**