

KIF3C

PDB:3B6V

Revision

Revision Type:created

Revised by:created

Revision Date:created

Entry Clone Accession:NP_002245

Entry Clone Source:

SGC Clone Accession:

Tag:N-terminal hexahistidine tag

Host:E.coli. BL21 (DE3) codon(+) RIL

Construct

Prelude:

Sequence:

mhahhhhhsgrenlyfqqASEALKVVARCRPLSRKEEAAGHEQILTMDVKGQVTLRNPRAAPGELPKTFTDAVYDASSKQADLYD
ETVRPLIDSVLQGFNGTVFAYGQTGTGKTYTMQGTWVEPELRGVIPNAFEHIFTTHISRSQNQQYLVRAASYLEIYQEEIRDLLSKEPG
KRLELKENPETGVYIKDLSSVTKNVKEIEHVMNLGNQTRAVGSTMNEVSSRSHAIFIITVECSERGSDGQDHIRVGKLNVLVDAG
SERQNKAGPNTAGGAATPSSGGGGGGGGGGAGGERPKEASKINLSALGNVIAALAGNRSTHIPYRDSKLTRLLQDSLGGNAKT
IMVATLGPASHSYDESLSTLRFANRAKNIKNKPRVNEDPKDTLLREF

Vector:p28a-mhl-tev

Growth

Medium:Terrific Broth

Antibiotics:

Procedure:We prepared the seeds by inoculating glycerol stock of E. coli cells BL21-CodonPlus (DE-3)-RIL containing the plasmid into 200 mL of Luria-Bertani medium. After overnight growth, all of the seeds were inoculated into 6 L of Terrific Broth medium in the presence of 50 μ g/mL of kanamycin and 50 μ g/mL chloramphenicol at 37 °C and grown to an OD600 between 3-5. Cells were then induced by isopropyl-1-thio-D-galactopyranoside at the final concentration of 0.5 mM and grown overnight at 18 °C in the SGC LEX bubbling system.

Purification

Procedure

1 ml of Ni-NTA suspension solution was added into 80 ml cell lysis supernatant solution. The mixture was shaken for 1 hour at 4 °C. Beads were collected with centrifuge at 2500 rpm, 5 minutes. Each fraction of Beads was washed with 100 ml washing buffer, then collected with centrifuge. Protein was eluted with 15 ml elution buffer.

The fractions eluted of the Ni-affinity chromatography applied to a Superdex S75 column

equilibrated in gel filtration buffer at a flow rate of 2.0 ml/min. Eluted fractions were 95% pure as judged by SDS-PAGE.

Extraction

Procedure

Cultures were centrifuged and the cell pellets were harvested and stored at -80 °C before use. Cells were thawed and suspended in 500 mL extraction buffer with 0.5% (v/v) protease inhibitor cocktail (Sigma), 1600 units Benzonase (Sigma), and lysed with microfluidizer. The lysate was centrifuged at 16000 rpm for 60 min and the supernatant was used for subsequent steps of purification. All the extraction steps were carried out at 4 °C.

Concentration: 28 mg/mL

Ligand

ADP and Mg(II)**MassSpec:**

Crystallization: Crystals were obtained using the vapor diffusion method and a protein concentration of 28 mg/ml containing 5 molar fold of ADP and MgCl₂. 2 μl of the concentrated protein mixed with 2 μl of a well solution containing 30% PEG 550 MME, 0.1 M (NH₄)₂SO₄, 0.1 M Na Cacodylate pH 6.5. Crystals appeared after two days at 18 °C.

NMR Spectroscopy:

Data Collection:

Data Processing: